

SUPPLEMENTARY TABLES

Please browse Full Text version to see the data of Supplementary Table 1

Supplementary Table 1. The results of the Proteome Profiler™ Array using the Mouse XL Cytokine Kit.

Supplementary Table 2. The primer sequences used in RT-PCR.

Fragment size (bp)	Name	Sequence (5' to 3')
342	HO-1-F	GCCCTGGAAGAGGAGATAGA
	HO-1-R	GTCGATGTTCCGGAAGGTAAA
288	SOD-1-F	CGGATGAAGAGAGGCATGTT
	SOD-1-R	GGTTTGAGGGTAGCAGATGAG
137	β-actin-F	ATCGTGCGAGACATCAATG
	β-actin-R	TCGTTGCCTATTGTGATGAC